






Contents

HOW TO ENTER 21CFR MODULE	3
HOW TO ADD USER	4
LOGIN USING ADMIN USER	5
HOW TO EDIT DEFAULT RIGHTS TO VARIOUS USER ROLES	6
HOW TO EDIT POLICY	8
HOW TO CHANGE PASSWORD	9
ADD SUPERVISOR / OPERATOR / QA USER	10



HOW TO ENTER 21CFR MODULE


- Start the SPAN Vision Studio Software. If it is operated for first time, login using
 - Username: spanuser Password: *****
- Administrator has to create his / her own user to login next time.
- Click 21CFR button as displayed in below image to enter 21CFR.





HOW TO ADD USER


- On entering 21CFR, click on "User" tab on right side of the screen.
- Click Add button on bottom side of the screen to add user.








- Enter the required information in "User Details" section.
- All the fields marked with asterisk are compulsory.
- Select User Role as Admin for first user.
- After all the information is entered, click "Generate" button to generate Electronic Signature.
- Click Save Button save newly created User.
- User details and rights can also be edited by click on.
- After Save system give a confirmation that user is created successfully.



Edit


LOGIN USING ADMIN USER

• Logout from current user "SPANUSER", by clicking on Logout Button as shown below.






• Re-login with new user which you have created.




HOW TO EDIT DEFAULT RIGHTS TO VARIOUS USER ROLES

- Again, go to 21CFR, click on Role tab.
- Click "Edit" button for role editing.





- Select Role to be edited.
- Select / Dis-select the desired rights to the selected role of user.
- Click save button once finished selecting rights.
- Now, whenever new user is created here after these new roles will be applied by default (however customization of rights is even possible during user management).




For example, select Alarm and Learn functionality to the Supervisor role which is dis-selected by default. (This is just for example)





HOW TO EDIT POLICY


- Click on Policy tab.
- Edit the required fields and click "Save".
- Note that this setting will be applicable to all the users.





HOW TO CHANGE PASSWORD

- Click Password tab.
- Enter Old Password, New Password and Confirm Password.
- Note that the new password complies with the password property set in security policy.





ADD SUPERVISOR / OPERATOR / QA USER

- Click User tab.
- Click "Add" button.
- Enter all the information in User Details.
- Select User Role as Supervisor / Operator / QA from drop down list.
- Click "Generate" button to generate Electronic Signature. This is allowed only to supervisor, QA and Admin users.
- You can select any of the previous users from "Configure from existing user", this will help you to assign rights to the selected user.

USER B = USER A (Can be easily done, you don't need to remember the rights given to USER A)

